



Development of diagnostic kits for selected markers of resistance, virulence and zoonotic transmission among methicillin-resistant *Staphylococcus aureus* strains

Dr. Boris Oberheitmann, Constanze Seidel
Q-Bioanalytic GmbH, BMBF-Projekt MetVetStaph II, IP11



Content

1. Background

- a) Q-Bioanalytic GmbH
- b) NALFIA

2. Developments

1. Triplex NALFIA test for *mecA*, *mecC* and *cyt b*
2. Improved DNA purification and optimization of sampling
3. Triplex CC398, *mecA* and *cyt b* Real-Time PCR

3. Outlook

Q-Bioanalytic GmbH

Existing portfolio



More than 20 Real-Time PCR kits
for food safety and quality testing



OneCup
Real-Time PCR



Selection of Real-Time PCR Kits

1. QuickBlue MRSA
2. OneCup *Salmonella*
3. OneCup *Listeria spec.*
4. QuickBlue *Listeria monocytogenes*
5. QuickBlue *Staphylococcus aureus*
6. QuickBlue *Cronobacter sakazakii*
7. QuickBlue *Clostridium perfringens*
8. QuickBlue *E. coli*
9. QuickBlue *Campylobacter jejuni*
10. QuickBlue *E. coli*, EHEC, EPEC, EIEC, *Shigella*
11. QuickBlue EHEC (stx1, stx2, eaeA)
12. QuickBlue *Vibrio vulnificus*
13. QuickBlue *Vibrio parahaemolyticus*
14. QuickBlue *Vibrio cholerae*
15. QuickBlue *Vibrio alginolyticus*
16. QuickBlue *Legionella pneumophila*
17. QuickBlue *Legionella spec.*
18. QuickBlue *Pseudomonas aeruginosa*

Q-Bioanalytic GmbH

Existing portfolio in medical microbiology



IVD CE marked
Real-Time PCR for MRSA



DNA Purification kit
based on magnetic
nano-particles






Q-Bioanalytic GmbH

Tasks:


- Development of PCR-based test kits using NALFIA and Real-Time PCR for rapid multianalyte diagnosis of resistance determinants
 - a) *mecA* and new *mecC* homologues
 - b) Differentiate between human and livestock-related MRSA lineages such as CC398
 - c) Typical relevant virulence or resistance markers

NALFIA


„Nucleic acid lateral flow immuno-assay (NALFIA)“ =
 Combining molecular biological principle of detection
 with immunochemical principle of visualization



NALFIA



Why NALFIA?




1. Bringing PCR to environments outside of labs and to Point-of-Care applications
2. Bringing the analytical knowledge to Point-of-Care users through smartphone applications

Developments

1. NALFIA test addressing *mecA*, *mecC* and *cyt b*

Detection of *mecA*, *mecC* gene and *cyt b* as internal amplification control (IAC) using NALFIA.

Detection limit: SCC*mec* VI (*mecA*) 1.5 pg, 10-100 cfu
SCC*mec* XI (*mecC*), 15 pg, 100-1000 cfu



DNA purification and optimization of sampling



Risk assessment concerning limit of detection in clinical samples revealed:

- Insufficient DNA purification can lead to false negative results
- Sampling with certain swabs can lead to insufficient release of the material into the purification reagents
- Time of analysis is a critical factor for acceptance of the methods in clinical settings


Conclusion:

- DNA extraction and/or sampling have to be optimized

Developments


Performance of prior existing DNA purification methods

QuickBlue DNA extraction kit (QBA):



LoD: 300 cfu

Extraction using spin column technology :



3×10^4 cfu

Developments

2. New optimized DNA purification procedure




Step 1
Swab sample
transferred into a tube
with lysis buffer
(thermal lysis)



Step 2
Supernatant
transferred to tube
with binding buffer
and silica-magnetite
nanoparticles



Step 3
Nano-particles are
magnetically immobilized
and re-suspended in 100 μ L
elution buffer




bioanalytic

Purification time max 30 min

Developments:

2. New DNA purification

Optimized
QuickBlue DNA
extraction kit
(QBA)



15000 cfu

1500 cfu

150 cfu

15 cfu

Optimization of sampling

Comparison of detection efficiency by Real-Time PCR

Classiq Swabs by COPAN:
Tip wrapped with
traditional Polyester
fiber.




Flocked Swabs by COPAN:
FLOQSwabs




Optimization of sampling,

Testing by dipping swabs in a decadal dilution series. A MRSA strain was grown 48hrs in Giolitti broth. Subsequently the culture was diluted and swabs were dipped into it in the presence and absence of 30 μ l blood.




Optimization of sampling

Classiq® Swabs by COPAN



FLOQSwabs® by COPAN



Results of experiments in the presence of blood



Results optimization of DNA purification and sampling

1. Risk of false negative results due to insufficient DNA purification could be reduced by reducing the limit of detection from 300 cfu to 15-150 cfu (30 minutes purification time)
2. Recovery rates from swabs in the presence of blood could be enhanced by an order of magnitude through application of FLOQSwabs




Developments

3. Triplex CC398, mecA and cyt b Real-Time PCR

- Specific primers for LA-MRSA CC398 after Stegger et al. (2011)
- Primers for the detection of mecA and cyt b
- Real-Time PCR with SybrGreen[®]
- Real-Time PCR with TaqMan[®] probe

Developments

3. Triplex CC398, *mecA* and *cyt b* Real-Time PCR




Developments

3. CC398 Real-Time PCR

Detection of *cc398* in
HEX channel

Detection of *mecA* in
Cy5 channel

Detection of *cyt b* in
FAM channel



3. Summary

Developments:

- a) Multiplex test for *mecA*, *mecC* gene and *cyt b* (IAC) using NALFIA
- b) Optimization of magnetic nano-particle-based purification method in combination with FLOQSwabs
- c) Triplex test cc398 with Real-Time PCR including internal control

Outlook:

Validation of the CC398 Real-Time PCR

4. Outlook

- Optimization work of the triplex Real-Time PCR system for detection of *cc398*, *mecA*, and *cyt b*
- Validation of the Real-Time PCR for detection of CC398:
 - Sensitivity
 - Specificity
 - Limit of detection
 - Robustness
- Development of NALFIAs to detect:
 - *mecA/C*, *S. aureus* and IAC
 - other MRSA resistance genes
 - MRSA virulence genes



Constanze Seidel
Dipl. Biol. Sonja Peters



PD Dr. rer. nat. Maren von
Köckritz-Blickwede



Prof. Dr. med. vet. Stefan Schwarz
Andrea T. Feßler, PhD

Vielen Dank für Ihre
Aufmerksamkeit!
Thanks for your co-operation